

## Sample Preparation Procedure for

# Fish and Seafood Trifluoroacetic Acid Extraction (TFAextr)

### General Information

This protocol for fish and seafood, provided by CVUA Karlsruhe, is based on the description of Mazzeo et al. [1] with modification, e.g. sample weight. It contains no purification step with magnetic beads. This method was used to create and validate a reference database for muscle meat of fish and seafood, which will be continuously expanded.

A selection of reference spectra (msp) and single spectra using this protocol is available for exchange on the MALDI-UP homepage (<https://maldi-up@ua-bw.de>).

### Field of Application

Muscle meat of fish, crustaceans and bivalve molluscs in natural form, raw, cooked or frozen.

### Chemicals and Material

- 1.5 ml reaction tubes
- Pipettes and tips for volumes from 1–500 µl
- Silica beads (0.1 mm resp. 0.5 mm diameter)
- Micropestle, compatible with the tubes used
- Bead mill (optional)
- Benchtop centrifuge
- HCCA matrix solution (please see »Tips and Recommendation«, page 2)
- Trifluoroacetic acid (TFA) 0.1 %

### References

- [1] Mazzeo, M. F., Giulio, B. D., Guerriero, G., Ciarcia, G., Malorni, A., Russo, G. L. & Siciliano R. A. (2008): Fish Authentication by MALDI-TOF Mass Spectrometry. *J Agric Food Chem* 56(23): 11071–11076.
- [2] Stephan, R., Jöhler, S., Oesterle, N., Näumann, G., Vogel, G. & Pflüger, V. (2014): Rapid and reliable species identification of scallops by MALDI-TOF mass spectrometry. *Food Control* 46: 6–9.

## Fish and Seafood: Trifluoroacetic Acid Extraction (TFAextr)

.....

### Extraction Procedure

- Add approx. 5 mg muscle without skin from the depth of the tissue into a 1.5 ml reaction tube
- Extract proteins either manually:
  - » Add 100 µl TFA (0.1 %) and approx. 100 mg silica beads into the tube
  - » Homogenize/grind the sample using a micropestle
  - » Add 400 µl TFA (0.1 %) and mixing thoroughly by vortexing for 10–20 sec. or mechanically:
  - » Add 500 µl TFA (0.1 %) into the tube
  - » Disrupt the cells using a bead mill with a frequency of 40 Hz for 3 minutes
- Leave for about 10 minutes (to solve proteins)
- Centrifuge at 15,900–18,400 rcf in a benchtop centrifuge for 10 minutes
- Pipet 1 µl supernatant onto a target sample spot (we recommend spotting the supernatant in duplicate or triplicate)
- As soon as the sample spot has dried, overlay the sample with 1 µl HCCA matrix solution (to prevent oxidation reactions which might cause peak shifts)
- Allow the sample spot to air dry before analysis → MALDI measurement

## Sample Preparation Procedure for

# Fish and Seafood Magnetic Beads Extraction (TFAextr C18)

### General Information

This protocol for fish, provided by CVUA Karlsruhe, is based on the description of Mazzeo et al. [1] with modification for e.g. sample weight and the purification step with magnetic beads. This method was used to create and validate a reference database for muscle meat of fish and seafood, which will be continuously expanded.

A selection of reference spectra (msp) and single spectra using this protocol is available for exchange on the MALDI-UP homepage (<https://maldi-up@ua-bw.de>).

### Field of Application

Muscle meat of fish, crustaceans, bivalve molluscs and fish eggs in natural form, raw, cooked or frozen.

### Chemicals and Material

- 1.5 ml reaction tubes
- Pipettes and tips for volumes from 1-500 µl
- Silica beads (0.1 mm resp. 0.5 mm diameter)
- Magnetic stand
- Micropestle, compatible with the tubes used
- Bead mill (optional)
- Benchtop centrifuge
- HCCA matrix solution (please see »Tipps and Recommendation«, page 2)
- Trifluoroacetic acid (TFA) 0.1 %
- Desorption solution (50 % acetonitrile in 0.1 % TFA)
- Magnetic beads (surface modified with C18 alkyl groups)

### References

- [1] Mazzeo, M. F., Giulio, B. D., Guerriero, G., Ciarcia, G., Malorni, A., Russo, G. L. & Siciliano R. A. (2008): Fish Authentication by MALDI-TOF Mass Spectrometry. *J Agric Food Chem* 56(23): 11071–11076.
- [2] Stephan, R., Jöhler, S., Oesterle, N., Näumann, G., Vogel, G. & Pflüger, V. (2014): Rapid and reliable species identification of scallops by MALDI-TOF mass spectrometry. *Food Control* 46: 6–9.
- [3] MoBiTec 2012, M-Beads-Proteomics S-C4, S-C8, S-C18

# Fish and Seafood:

## Magnetic Beads Extraction (TFAextr C18)

.....

### Protein Extraction

#### 1. Protein Extration

- Add approx. 5 mg muscle without skin from the depth of the tissue into a 1.5 ml reaction tube
- Extract proteins either manually:
  - » Add 100  $\mu$ l TFA (0.1 %) and approx. 100 mg silica beads into the tube
  - » Homogenize/grind the sample using a micropestle
  - » Add 400  $\mu$ l TFA (0.1 %) and mixing thoroughly by vortexing for 10–20 sec. or mechanically:
  - » Add 500  $\mu$ l TFA (0.1 %) into the tube
  - » Disrupt the cells using a bead mill with a frequency of 40 Hz for 3 minutes
- Leave for about 10 minutes (to solve proteins)
- Centrifuge at 15,900–18,400 rcf in a benchtop centrifuge for 10 minutes

### Purification of Protein Extract with Magnetic Beads

#### 2. Conditioning / Resuspension Magnetic Beads

- Pipette 20  $\mu$ l of the magnetic beads solution into a 1.5 ml reaction tube, place the tube on the magnet, remove the supernatant with a pipette and discard it
- Remove the reaction tube from the magnet, add 100  $\mu$ l TFA (0.1 %) and resuspend the magnetic beads
- Place the tube on the magnet, remove the supernatant with a pipette and discard it, repeat this washing step a second and a third time
- Add 10  $\mu$ l TFA 0.1 % and resuspend the magnetic beads

#### 3. Protein/Peptide adsorption

- Add 25  $\mu$ l sample solution from step 1 into the tube containing magnetic beads and leave for 10 minutes at ambient temperature
- Place the tube on the magnet, remove the supernatant with a pipette and discard it
- Remove the reaction tube from the magnet, add 50  $\mu$ l TFA (0.1 %) and resuspend the magnetic beads
- Place the tube on the magnet, remove the supernatant with a pipette and discard it, repeat this washing step a second and a third time

#### 4. Protein/Peptide desorption

- Add 10  $\mu$ l desorption solution to the sample-coated magnetic beads from step 3, resuspend and incubate for 8 minutes at ambient temperature
- Pipet 1  $\mu$ l supernatant onto a target sample spot (we recommend spotting the supernatant in duplicate or triplicate)
- As soon as the sample spot has dried, overlay the sample with 1  $\mu$ l HCCA matrix solution (to prevent oxidation reactions which might cause peak shifts)
- Allow the sample spot to air dry before analysis → MALDI measurement