

Sample Preparation Procedure for

Diary Products Organic Solvent Extraction (OSExtr)– (incl. Variant with Protein Precipitation)

General Information

This protocol for dairy products was provided by CVUA Stuttgart and used to create and validate a reference database for feta and mozzarella cheese [1, 2].

A comprehensive collection of reference spectra (msp) and single spectra generated using this protocol is available for exchange on the MALDI-UP homepage (<https://maldi-up@ua-bw.de>).

Field of Application

Dairy products, e.g. milk, yogurt or cheese.

For pure milk it is necessary to perform protein precipitation (steps 1–4).

For yoghurt or (fresh-)cheese (already “precipitated” casein) please start with step 5.

Chemicals and Material

- 2.5 ml reaction tubes and tips.
- Spatula
- Centrifuge with centrifuge tubes and temperature control
- Small beaker
- Heating block
- Benchtop centrifuge
- Dilute acetic acid (50%)
- Aqua dest/bidest
- HCCA Matrix solution (please see »Tipps and Recommendation«, page 2)
- OS Solvent (acetonitrile 50%, water 47.5% and trifluoroacetic acid 2.5%)

References

- [1] Männig, A., Hiller, E., Rau, J. (2017): MALDI-TOF MS zur Tierartenbestimmung bei Milch und Käse. 46; Poster at the “Deutscher Lebensmittelchemikertag” 25.–27.09.2017, Würzburg. https://maldi-up.ua-bw.de/docs/CVUAS_Maennig_et_al_Tierarten_MALDITOF_2017.pdf
- [2] Rau, J., Korte, N., Dyk, M., Wenninger, O., Schreiter, P., Hiller, E. (2020a): Rapid Animal Species Identification of feta and mozzarella cheese using MALDI-TOF mass-spectrometry. Food Control, 107349. <https://doi.org/10.1016/j.foodcont.2020.107349>

Diary Products: Organic Solvent Extraction (OSExtr) – incl. Variant with Protein Precipitation

.....

Protein / Peptides Precipitation

1. Separation of Fat/Skimming

For milk with especially low fat content, e.g. mare milk, this step can be skipped.

- Add approx. 30 ml milk to a centrifuge tube
- centrifuge at 3,345 rpm and 8°C for 12 minutes
- To remove fat, please carefully decant the liquid into a small beaker; the fat should cling to the wall of the centrifuge tube and is to be discarded.

2. Precipitation of Protein / Peptides (mostly casein)

- Add a few drops of dilute acetic acid to the skimmed milk
- Heat the liquid to 40–45°C (e.g. in a heating block) for at least 15 minutes until sufficient casein is precipitated.

3. Separation of Casein and Whey

- Transfer the liquid with the precipitate to 1.5 ml reaction tubes

For mare milk you will need to pool the content of about 6 Eppendorf tubes, since mare milk does not seem to contain as much casein as other milk types.

For other milk types less tubes are required

- Centrifuge at 12,000–14,000 rcf in a benchtop centrifuge for 2 minutes
- Discard the supernatant (= whey) and merge the pellets (=casein)

4. Washing

- Add 1 ml aqua dest. and stir with a spatula
- Centrifuge at 12,000–14,000 rcf in a benchtop centrifuge for 2 minutes
- Discard the supernatant

Extraction Procedure

5. Organic Solvent Extraction Procedure

- Transfer a small amount (about 3 mm³) of the sample or protein precipitate into a 1.5 ml tube
- Add 200 µl OS
- Homogenize by stirring with a spatula
- Mix thoroughly by vortexing for 10–20 seconds
- Centrifuge at 12,000–14,000 rcf in a benchtop centrifuge for 2 minutes
- Pipet 1 µL supernatant onto a target sample spot (We recommend spotting the supernatant in duplicate or triplicate)
- As soon as the sample spot has dried, overlay the sample with 1 µL HCCA matrix solution (to prevent oxidation reactions leading to peak shifting)
- Allow the sample spot to air dry before analysis → MALDI measurement