

MSP creation for MALDI Biotyper



| |
|---|
| ID (culture collection number or similar) |
|---|

Metadata (MSP-Metadata MBT Compass):

| | | | | | |
|---|--------------------------|--|-------------------------------|--------------------------------|---------------------|
| Organism | | | | | |
| Strain (e.g. ATCC nr./ ID/ ...) | | | | | |
| Provided by (e.g. ATCC/DSMZ/ ...) | | | | | |
| Determined by (sequenced/ type strain/ ...) | | | | | |
| Conserved | <input type="checkbox"/> | Sample Preparation ("Extraction Method") | <input type="checkbox"/> DT | <input type="checkbox"/> eDT | |
| | | | <input type="checkbox"/> EFex | <input type="checkbox"/> | |
| Matrix | | HCCA | | | |
| Growing conditions | | Agar | Temperature (°C) | Time (h) | Culture Conditions: |
| Comment | | | | | |

Information of measured sample data:

| | | | | |
|---------------------------------------|--------------------------------------|-------|-------|----------|
| Count of measured spectra | | Date: | Time: | Acronym: |
| File directory of the raw data | <input type="checkbox"/> D/Data /... | | | |

Spectra editing (flexAnalysis):

- Load the measured sample spectra and the BTS
- **WINDOWS EXPLORER: rename ...**
 - file „BTS“ → „BTS raw“
 - raw sample spectra file: e.g.: ID 1234 → ID 1234 **raw 24sp**

| | | | |
|---|--------------------------|--------------------------|--------------------------|
| Select all spectra → | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| Assign Method Method → Open ... | MBT_Standard.FAMSMETHOD | Baseline Subtraction | Smooth (1x) |

BTS check/ recalibration:

| | | | |
|--|--|--|--|
| Check Mass Control List Calibrate → Internal ... | <input type="checkbox"/> Automatic-Assign | | <input type="checkbox"/> Peaks manually assigned |
| | Max. deviation (ppm): | | |
| Recalibrate sample spectra | <input type="checkbox"/> Copy calibration | | |
| Mass calibration constants BTS Select BTS spectrum → Properties ... | C0: | | |
| | C1: | | |
| | C2: | | |
| Mass calibration constants sample spectra | <input type="checkbox"/> Check: same as BTS? | | |

- Close and save the BTS spectrum
 - WINDOWS EXPLORER: rename** the newly created file (by flexAnalysis) „BTS“ → „BTS ed“

Editing the sample spectra:

| | | | | | | |
|---|--------|--------|--------|--------|---------|----------------------------------|
| Conspicuous spectra (position\measurement): (Flat lines etc.) | | | | | | <input type="checkbox"/> removed |
| Remaining spectra: Peak accuracy (calculation Excel-worksheet, +/- 500ppm) | | | | | | |
| m/z | ≈ 3000 | ≈ 5000 | ≈ 6000 | ≈ 8000 | ≈ 10000 | |
| Minimum Mass (top of the peak(s)) | | | | | | |
| Maximum Mass (top of the peak(s)) | | | | | | |
| Removed spectra | | | | | | |
| Count of remaining spectra | | | | | | |

- **Select removed spectra and close** (right click → „Close“) → **DO NOT SAVE THE SPECTRA!**
- **Close remaining spectra and SAVE THEM ALL!**

| | |
|--|-------------------|
| WINDOWS EXPLORER: rename the newly created file (by flexAnalysis): e.g. ID 1234 → ID 1234 <i>ed 21sp</i> | File name: |
|--|-------------------|

MSP Creation with MBT Compass Explorer:

- Open MBT Compass Explorer
- Load and select all edited sample spectra
- Right click → „Create MSP“ → **assign MSP name**

| |
|--|
| MSP Name: e.g.: Streptobacillus hongkongensis DSM 26322 CVUAS / Escherichia coli CVUAS 5146 CVUAS |
|--|

- **Taxonomy tree:** change the dropdown list to “Projects”, select the node where the MSP should be stored and start the Taxonomy Tree Editor (right click or button next to dropdown menu)

Edited MSP **Metadata**

Added MSP to “Projects” node:

Verification of the MSP with an **independent spectrum** (date):

Report printout

Preparation:

- DT
- eDT
- EFex
-

MSP created ...: own MSP library updated

| | |
|-----------------|--|
| Comment: | |
|-----------------|--|

Date / acronym _____